

*As we enter our fifth year, iRad Technical Services is happy to introduce Kodak's new DirectView CR Mammography Feature and Kodak's entry into clinic DR market with the DirectView DR 3000.*

## **Introducing CR Mammography from Kodak**

If you already use the CR 850 or the CR 950/975, this feature allows you to use the system in dedicated mammography departments as well as mixed-use environments.

**It requires no change in CR procedure, workflow or user interface and complements the CR investment that you have already made.**

### **Mammography-quality images**

Kodak builds on its impressive track record for delivering exceptional image quality with the Kodak DirectView CR system. See high-resolution scans with a 50-micron pixel pitch and more consistent image quality than conventional imaging.

Kodak DirectView CR Mammo cassettes allow for image capture close to the chest wall. The new EHR-M2 screen delivers the sharpest CR mammography screen.

In addition, Kodak's DirectView EVP Software and Kodak Black Surround/Masking Software provide optimal quality for mammography.

The CR Mammography Feature is an option for purchase with a new CR 850 or CR 950/975, or as an upgrade to an existing system.

The single-cassette CR 850 is ideal for distributed applications. The CR 975 system manages 16 cassettes simultaneously, providing centralized processing with "drop-and-go" decentralized workflow.


### **CR Mammography Total Quality Tool**

*The new Mammography Feature even provides an efficient and accurate way to verify and report mammography image quality!*

The Kodak DirectView CR Mammography Total Quality Tool is a complete package that includes everything you need for easy, accurate and objective mammography image testing. It enables you to perform image tests and QC measurements with the same interface used for examinations.


Tests are done at your facility and at your convenience — without the cost, scheduling or disruption of third-party testing.

For more information, call **iRad at 905.844.3897 or 1-877-844-3810**, or go online at [www.kodak.com/go/mammo](http://www.kodak.com/go/mammo)



## Introducing Kodak DirectView DR Systems

Kodak Digital Radiography systems utilize a high-performance image capture chain for converting X-rays into digital images.



### *DirectView DR 3000 for General Radiography:*

*Kodak's DirectView DR 3000 is a powerful multipurpose system designed to handle a full range of general radiography procedures.*

It is capable of handling a variety of image capture options including Upright, Table and Angled.

In this fully motorized system, the arm maintains constant alignment between the x-ray tube and the detector for easy and fast patient positioning. Some of the features incorporated into the DR 3000 are:

- User-friendly control panels which are conveniently located at both the Bucky and the tube side of the arm.
- Grids that are easily removed and a large image area, 43 x 43 cm, which eliminates the need to rotate the detector.
- Collimator filters that can be selected with the push of a button and emergency system stop buttons and sensors that help assure safety of operators and patients.

### *The space-saving design of the DR 3000 System ensures that it can be placed in almost any room*

The floor-mounted stand holds both the x-ray tube and the Bucky. This freestanding system requires no ceiling suspension or wall mounts.

### *Superb image quality with preview images available in just 8 seconds!*

For more information, call **iRad** at 905.844.3897 or 1-877-844-3810, or go online at [www.kodak.com/go/dr](http://www.kodak.com/go/dr)

**iRad Technical Services Inc.**  
1235 Trafalgar Road, Suite 306  
Oakville, Ontario  
L6H 3P1

Tel: 905.844.3897 Toll-free: 877.844.3810  
Fax: 905.845.7161